## IK-161 (F/A 750581)

INSTALLATION INSTRUCTIONS FOR FILTRATION ADAPTER KIT USED ON GLASS TUBES WITH TARGET ANGLES LESS THAN 16° IN B-100, B-130, OR B-150 HOUSINGS

## IK-162 (F/A 750580)

INSTALLATION INSTRUCTIONS FOR FILTRATION ADAPTER KIT USED ON GLASS TUBES WITH 16° TARGET ANGLES IN B-100, B-130, OR B-150 HOUSINGS

### DESCRIPTION

The IK-161 or IK-162 is an adapter kit designed to allow the installation of multiple aluminum discs of varying thickness such that the beam filtration without the collimator can be adjusted from insert inherent filtration level to 3.0mm of Aluminum maximum. (Value of filter discs: 0.5mm, 0.8mm, & 1.0mm)

The intent of this adapter kit is to initially provide a total filtration of 1.5mm Aluminum at the tube unit (0.7mm inherent plus 0.8mm added) along with a label indicating 1.5mm total filtration. The installation company is then responsible for the labeling and installation of other filter values. The installation company is also responsible for local regulations regarding total filtration.

### INSTALLATION TO MEET 1.5mm TOTAL FILTRATION:

CAUTION - DO NOT REMOVE PORT RETAINING PLATE!

USE PROPER RETAINING RING TOOLS AND EXERCISE EXTREME CARE TO NOT DAMAGE THE Pb RETAINING SHOULDER.

- Refer to page 2 drawing for part location
- Remove item (1) Port Cover Plate
- Remove Pb Beam Diaphragm
- Install item (2) Pb Filter Cone
- Replace item (1) Port Cover Plate
- Install the 0.8mm Aluminum Disc
- Install item (4) Wave Spring
- Install item (5) Retaining Ring
- Apply 1.5mm Total Filtration label directly under the tube unit Minimum Filtration label.

### INSTALLATION TO PROVIDE ADDITIONAL FILTRATION BEYOND THE INITIAL 1.5mm TOTAL:

CAUTION - DO NOT REMOVE PORT RETAINING PLATE!

USE PROPER RETAINING RING TOOLS AND EXERCISE EXTREME CARE TO NOT DAMAGE THE Pb RETAINING SHOULDER.

- Remove item (5) Retaining Ring
- Remove item (4) Wave Spring
- Install Aluminum discs to achieve total filtration requirement (0.5mm and 1.0mm available)
- Install item (4) Wave Spring
- Install item (5) Retaining Ring
- Apply Total Filtration label provided by installer

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