Note: Document originally drafted in the English language.  
注释: 文件最初用英语起草。

### Product Description

The MCS-8064 is a 9.4" (240 mm) 140 kV, 5.7 MJ (8.0 MHU) maximum anode heat content, rotating anode insert. This insert is specifically designed for GE CT Scanners. The insert features a 7° tungsten-rhenium facing on molybdenum with a graphite backed target and is available with the following nominal focal spots values:

0.9 x 0.7  
1.2 x 1.2

**Maximum Anode Cooling Rate:**  
21,800 W (30,520 HU/sec)

**Maximum continuous anode heat dissipation:**  
8,000 W (11,200 HU/sec)

**Nominal Anode Input Power:  
Load Time 1 Second**  
Small - 47 kW IEC 60613:2010  
Large - 100 kW IEC 60613:2010

**Nominal CT Anode Input Power:**  
Small - 47 kW IEC 60613:2010  
Large - 96 kW IEC 60613:2010

**Reference Axis:**  
Perpendicular to port face.

This insert is intended for use in the Varex Imaging B-680H housing.

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The single exposures are controlled by system software.

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† Available only with VCT 85 kW Option installed. Otherwise constrained to the max mA available.

\* Available only with VCT 85 kW & VCT Hi Power Options installed. Otherwise constrained to the max mA available.

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The serial exposures are controlled by system software and are applicable for repeat every 10 minutes based on 3 hour wait period after tube warm up has completed.

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† Available only with VCT 85 kW Option installed. Otherwise constrained to the maximum mA available.  
 \* Available only with VCT 85 kW & VCT Hi Power Options installed. Otherwise constrained to the maximum mA available.

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The system can acquire 72 kW at 120 kVp for 4 seconds scan duration. The system can acquire 84 kW at 120 kVp for 4 seconds scan duration if the 85 kW Option is enabled. The system can acquire 96 kW at 120 kVp for 4 seconds scan duration if the 85 kW and Hi Power Options are enabled. The single exposures are controlled by system software.

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



1.2 x 1.2 

**Product Description**

Maximum Peak Voltage ..... 140 kV  
 Anode to Ground ..... 0 kV  
 Cathode to Ground ..... 140 kV

Maximum X-ray Tube Assembly Heat Content  
 ..... 5.5 MJ (7.7 MHU)

Maximum Housing Temperature ..... 78°C

Maximum Continuous Heat Dissipation @ 25°C ambient  
 (Includes stator heat) ..... 8.0 kW (11.2 kHU/sec)

Maximum Heat Exchanger Dissipation  
 ..... 8.0 kW (11.2 kHU/sec)

Permanent Filtration  
 X-Ray Tube Assembly (IEC 60522) 3.25mm Al / 70kV  
 X-Ray Insert ..... .02mm Al / 70kV

Loading Factors for Leakage Radiation  
 ..... 140 kV, 57 mA

High Voltage Cable ..... Special

Ambient Air Temperature Limits for Operation  
 ..... 15°C to 45°C

Temperature Limits for Storage and Transport  
 ..... -20°C to + 75°C  
 Humidity ..... 10% to 90%  
 Atmospheric Pressure Range ..... 70 kPa to 106 kPa

Weight - Housing ..... 106 kg (234 lbs)

IEC Classification ..... Class 1

Safety Devices  
 Housing - Thermal Switch  
 Normally Closed Contacts ... Opens at 93°C ±3°C

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




Dimensions are for reference only  
维度是供仅参考

## Tube Housing Assembly Heating and Cooling 管套组件加热与冷却


**Note:**

1. Heat inputs into housing include tube power, filament power, and stator power.
2. Heating curves based on no restrictions of natural convection around tube housing assembly.
3. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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Terminal / Wire Color Chart  
终端 / 导线颜色图


<b>Stator Type:</b> 3 Ø, 4 Pole	
<b>Stator Coil Resistance:</b>	2.3 Ohms ± 5%
<b>Starter Voltage:</b>	<b>Start</b> <b>Run</b>
	280 Hz      430 VAC    260 VAC
<b>Time to Full Speed:</b>	
140 Hz ±2 Hz in < 12 Sec.	
8400 RPM ±1120 RPM	
<b>X-Ray Tube Assembly:</b>	
MCS-8064/B-680H	IEC 60601-2-28

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Anode Heating and Cooling Curves  
阳极加热与冷却曲线

**Note:**  
Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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