Note: Document originally drafted in the English language.
Note : Document à l'origine rédigé dans l'anglais.
Anmerkung: Dokument ursprünglich gezeichnet in der englischen Sprache.
Nota: Documento elaborado originalmente en la lengua inglesa.

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The single exposures are controlled by system software.

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NOTE: 72 kW base configuration is not available in all markets.

Large Spot Single Exposure Limits				
Scan Time	140 kV	120 kV	100 kV	80 kV
5	615*	800*	760*	675 [†]
10	575 [†]	790*	760*	675 [†]
20	535 [†]	645 [†]	760*	675 [†]
30	490	560	695 [†]	675 [†]
40	445	510	635 [†]	675 [†]
50	410	485	585	675 [†]
60	385	460	550	675 [†]

[†] Available only with VCT 85 kW Option installed. Otherwise constrained to the max mA available.
^{*} Available only with VCT 85 kW & VCT Hi Power Options installed. Otherwise constrained to the max mA available.

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The serial exposures are controlled by system software and are applicable for repeat every 10 minutes based on 3 hour wait period after tube warm up has completed.

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† Available only with VCT 85 kW Option installed. Otherwise constrained to the maximum mA available.
 * Available only with VCT 85 kW & VCT Hi Power Options installed. Otherwise constrained to the maximum mA available.

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The system can acquire 72 kW at 120 kVp for 4 seconds scan duration. The system can acquire 84 kW at 120 kVp for 4 seconds scan duration if the 85 kW Option is enabled. The system can acquire 96 kW at 120 kVp for 4 seconds scan duration if the 85 kW and Hi Power Options are enabled. The single exposures are controlled by system software.

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

THREE PHASE EMISSION (± .15 A)
 0.9 x 0.7 

THREE PHASE EMISSION (± .15 A)
 1.2 x 1.2 





Maximum Peak Voltage	140 kV
Anode to Ground	0 kV
Cathode to Ground	140 kV
Maximum X-ray Tube Assembly Heat Content	5.5 MJ (7.7 MHU)
Nominal Continuous Input Power @ 25°C ambient	8.0 kW (11.2 kHU/sec) IEC 60613:2010
Maximum Housing Temperature	78°C
Maximum Heat Exchanger Dissipation	8.0 kW (11.2 kHU/sec)
Permanent Filtration	
X-Ray Tube Assembly (IEC 60522)	3.25mm Al / 70kV
X-Ray Insert	.02mm Al / 70kV
Loading Factors for Leakage Radiation	140 kV, 57 mA
High Voltage Cable	Special
Ambient Air Temperature Limits for Operation	15°C to 45°C
Temperature Limits for Storage and Transport	-20°C to + 75°C
Humidity	10% to 90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight - Housing	106 kg (234 lbs)
IEC Classification	Class 1
Safety Devices	
Housing - Thermal Switch	
Normally Closed Contacts	Opens at 93°C ±3°C
Filament Frequency Limits	50 Hz - 25 kHz
Power Supply	DC

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
MCS-8064
TUBE HOUSING ASSEMBLY HEATING AND COOLING

Note:

- Heat input into housing includes all power sources; tube, filament, stator and circulating pump.
- Heating curves based on no restrictions to air flow through heat exchanger, or natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

Note:
 Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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