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
Small - White  
小焦点 - 白色

Stand - By  
备用

Frame or Chasis  
框架或底盘

X-Ray Tube  
X 射线管

Radiation Filter or Filtration  
辐射过滤器或过滤

Dimensions are for Reference only  
维度是供仅供参考

Note: Document originally drafted in the English language.

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### Product Description

The RAD-94 is a 4"(100 mm), 150 kV, 450 kJ (600 KHU) rotating anode insert specifically designed for general radiography, fluoro/spotfilm radiography, simulator work for CT extension applications and when using fluoroscopic digital systems. The insert features a 14° tungsten rhenium molybdenum graphite target and is available in the following focal spot combinations:

0.4 - 0.8  
IEC 60336

### Nominal Anode Input Power

Small - 15 kW IEC 60613  
Large - 45 kW IEC 60613

For the equivalent anode input power of 190 Watts

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



50 Hz - 2,850 RPM

60 Hz - 3,450 RPM

For 1Ø and other applications, please consult the manufacturer

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



## CINERADIOGRAPHIC RATINGS HOW TO USE CINERADIOGRAPHIC CHARTS

**General:** With the Cineradiographic rating chart we can determine the maximum allowable kW of the Cine pulse, or with a given kW determine maximum time in seconds the Cine run can progress.

The Most common way of using the charts is to determine maximum time of any expected Cine run and maximum duty factor. With a known duty factor and Cine run time kW can easily be determined.

### Definition of Terms

**Time in seconds:** Total time of one Cine run, usually 5 to 12 seconds.

**Duty Factor in Percent (DF%):** Actual time during one second the x-ray tube is producing x-rays. If we select a 4 msec pulse width and 60 exposures per second the x-ray tube will be producing x-rays for a total of 240 msec each second or 24% of the time. The higher the DF number, the more load placed on the x-ray tube.

**Peak Pulse Power:** Peak energy in watts of any one Cine Pulse. Can be any combination of kV and mA allowed by Radiographic and Filament Emission curves.

Example: 80 kV at 400 mA equals

$$80,000 \text{ V} \times 0.4 \text{ A} = 32,000 \text{ W or } 32 \text{ kW}$$

### USING THE CINE RATING CHARTS:

RAD-94 150/180 Hz 3 Phase 0.8 Focal Spot

**Example:** Determine maximum kW allowed with the following known factors:  
 Maximum Pulse Width ..... 4 msec  
 Exposures per Second ..... 60  
 Maximum Cine Run Time ..... 10 seconds

### Calculate Duty Factor: (DF%)

$$DF\% = \frac{\text{Pulse Width (mSec)} \times \text{Frames per Second}}{10}$$

$$DF\% = \frac{4 \text{ msec} \times 60 \text{ exp/sec}}{10} = \frac{240}{10} = 24\%$$

Refer to Rating Chart RAD-94 150/180 Hz 3 Phase 0.8 Focal Spot:

At bottom of chart find 10 second line. Move vertically to intersection with 24% DF curve. Make a horizontal reference to left side of rating chart and note kW rating of 34 kW.

We now know each pulse during the cine run can have a maximum rating of 34 kW under conditions given in example.

kW = kV x mA. The kW of the exposure can be any combination of mA and kV allowed by the Radiographic and Filament Emission Charts.

The Cine rating charts are usable to 100% anode heat storage. The start of Cine run should be below 70% and heat storage. Exceeding 100% anode heat storage will cause anode track erosion with high risk of tube destruction.

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Manufactured by Varex Imaging Corporation